



PCRBIOSYSTEMS

simplifying research

PCR BIO Classic Taq

Product description:

PCR BIO Classic Taq is a highly-purified, recombinant Taq DNA Polymerase. The latest developments in buffer chemistry allow for enhanced PCR speed, yield and specificity. The enzyme and buffer system allow for superior PCR performance on complex templates such as mammalian genomic DNA.

PCR BIO Classic Taq is a robust enzyme for all your everyday PCR applications including genotyping, screening and library construction. PCR BIO Taq DNA Polymerase performs consistently well on a broad range of templates (including both GC and AT rich). PCR BIO Classic Taq has 5'-3' exonuclease activities, but no 3'-5' exonuclease (proofreading) activity.

The enzyme has an error rate of approximately 1 error per 2.0×10^5 nucleotides incorporated. PCR products generated with PCR BIO Classic Taq are A-tailed and may be cloned into TA cloning vectors.

Component	1500 units	3000 units	6000 units
PCR BIO Classic Taq (5u/ μ l)	1x 300 μ l	2x 300 μ l	4x 300 μ l
10x PCR BIO Classic Buffer	3x 2ml	6x 2ml	12x 2ml
50mM MgCl ₂	2x 2ml	4x 2ml	8x 2ml

Shipping and Storage

On arrival the kit should be stored at -20°C. Avoid prolonged exposure to light. If stored correctly the kit will retain full activity for 12 months. The kit can be stored at 4°C for 1 month. The kit can go through 30 freeze/thaw cycles with no loss of activity.

Limitations of product use

The product may be used only for in vitro research purposes.

Technical support

For technical support and troubleshooting please email technical@pcrbio.com the following information:

- Amplicon size
- Reaction setup
- Cycling conditions
- Screen grabs of gel images

www.pcrbio.com

Important considerations

PCRBIO 10x reaction buffer: The 10x reaction buffer contains enhancers and stabilizers. It is not recommended to add further PCR enhancers to the reaction. The buffer composition has been optimised to maximise PCR success rates.

MgCl₂ and dNTPs: A final reaction concentration of 3mM MgCl₂ to 1mM dNTPs (0.25mM each) is recommended.

Template: For eukaryotic DNA use between 5ng and 500ng per reaction, for cDNA use below 100ng per reaction.

Primers: Primers should have a predicted melting temperature of around 60°C, using default Primer 3 settings (<http://frodo.wi.mit.edu/primer3/>). The final primer concentration in the reaction should be between 0.2µM and 0.6µM.

Annealing: We recommend performing a temperature gradient to experimentally determine the optimal annealing temperature. Alternatively, we recommend a 55°C annealing temperature then increase in 2°C increments if non-specific products are present.

Extension: Optimal extension is achieved at 72°C. The optimal extension time is dependent on amplicon length and complexity of template. 15 seconds per kilobase (kb) is recommended for amplification from eukaryotic DNA for amplicons between 1kb and 6kb. For shorter amplicons a 1 second extension is sufficient.

Reaction setup

1. Prepare a master mix based on the following table:

Reagent	50µl reaction	Final concentration	Notes
10x PCRBIO Classic Buffer	5.0µl	1x	
100mM dNTPs (25mM each)	0.5µl	1mM (0.25mM each)	
50mM MgCl ₂	3.0µl	3mM	
Forward primer (10µM)	2.0µl	400nM	See above for optimal primer design
Reverse primer (10µM)	2.0µl	400nM	See above for optimal primer design
Template DNA	<100ng cDNA, <500ng genomic	variable	See above for template considerations
PCRBIO Classic Taq (5u/µl)	0.25µl - 1.0µl		

2. Cycle using conditions based on the following table:

Cycles	Temperature	Time	Notes
1	95°C	1min	Initial denaturation
40	95°C	15 seconds	Denaturation
	55°C to 65°C	15 seconds	Anneal
	72°C	1 to 90 seconds	Extension (15 seconds per kb)