



**PCRBIOSYSTEMS**

simplifying research

## qPCR BIO cDNA Synthesis Kit

### Product description:

The qPCR BIO cDNA Synthesis Kit uses the latest developments in reverse transcriptase technology and buffer chemistry to enhance cDNA synthesis speed and yield with accurate transcript representation. The reverse transcriptase, buffer system and combination of random hexamers with anchored oligo(dT) allow for unbiased, efficient, sensitive cDNA synthesis.

The modified MMLV reverse transcriptase (RTase) is both thermostable and extremely active. The RTase is not inhibited by ribosomal and transfer RNAs, total RNA is an ideal substrate. The enzyme is blended with RNase inhibitor preventing degradation of RNA by contaminating RNase.

5x buffer contains anchored oligo(dT), random hexamers, enhancers, dNTPs and MgCl<sub>2</sub>. The relative concentrations of random hexamers and anchored oligo(dT) have been optimised for the generation of cDNA for use in real-time PCR experiments. The kit can be used with 4.0pg to 4.0µg total RNA.

| Component                        | 25 reactions | 100 reactions | 500 reactions |
|----------------------------------|--------------|---------------|---------------|
| 5x cDNA synthesis mix            | 100µl        | 400µl         | 2000 µl       |
| 20x RTase (with RNase inhibitor) | 25µl         | 100µl         | 500µl         |

### Shipping and Storage

On arrival the kit should be stored at -20°C. Avoid prolonged exposure to light. If stored correctly the kit will retain full activity for 12 months. The kit can be stored at 4°C for 1 month. The kit can go through 30 freeze/thaw cycles with no loss of activity.

### Limitations of product use

The product may be used only for in vitro research purposes.

### Technical support

For technical support and troubleshooting please email [technical@pcrbio.com](mailto:technical@pcrbio.com) the following information:

Reaction setup  
PCR cycling conditions  
Screen grabs of gel images / real-time PCR traces

[www.pcrbio.com](http://www.pcrbio.com)

## Important considerations

**5x cDNA Synthesis Mix:** Contains anchored oligo(dT), random hexamers, 15mM MgCl<sub>2</sub>, 5mM dNTPs, enhancers and stabilizers. It is not recommended to add further enhancers or MgCl<sub>2</sub> to the reaction. The buffer composition has been optimised to generate cDNA for downstream real-time PCR analysis.

**Template:** For total RNA use between 4pg and 4.0µg per reaction.

**Incubation temperature:** We recommend incubating with a temperature of 42°C for 30 minutes for the majority of applications (<65% GC). Where regions of interest contain high secondary structure (>65% GC) incubation temperatures of up to 55°C may be used.

**qPCR setup:** We recommend 4.0µl of cDNA per 20µl real-time PCR reaction.

## Reaction setup

1. Allow 5x cDNA Synthesis Mix to thaw, briefly vortex.
2. Prepare a master mix based on the following table. Insert reagents in sequence listed:

| Reagent                           | 20µl reaction           | Final concentration | Notes   |
|-----------------------------------|-------------------------|---------------------|---|
| 5x cDNA synthesis mix             | 4.0µl                   | 1x                  |   |
| 20x RTase                         | 1.0µl                   |                     | Add before total RNA as RNase inhibitor is blended with RTase |
| Total RNA (between 4pg and 4.0µg) | Xµl                     |                     | Variable  |
| PCR grade dH <sub>2</sub> O       | Up to 20µl final volume |                     |   |

## No RT control setup (optional)

3. Prepare a master mix based on the following table. Insert reagents in sequence listed:

| Reagent                           | 20µl reaction           | Final concentration | Notes                         |
|-----------------------------------|-------------------------|---------------------|-------------------------------|
| 5x cDNA synthesis mix             | 4.0µl                   | 1x                  |                               |
| Total RNA (between 4pg and 4.0µg) | Xµl                     |                     | Use equal amount as in step 2 |
| PCR grade dH <sub>2</sub> O       | Up to 20µl final volume |                     |                               |

## Incubation and enzyme denaturation

4. Incubate at 42°C for 30 minutes.
5. Incubate at 85°C for 10 minutes to denature RTase.