



PCRBIOSYSTEMS

simplifying research

PCR BIO HS Taq 500units/ μ l

Product description:

PCR BIO Taq (500u/ μ l) is a robust enzyme for all your everyday PCR applications including genotyping, screening and diagnostics. Glycerol, detergents and many PCR enhancers are often prohibitive for successful lyophilisation. The use of a highly concentrated enzyme preparation ensures that components prohibitive to lyophilisation are diluted out. The hot start (HS) additive is supplied in 10mM Tris pH 8.0, it does not contain glycerol.

PCR BIO Taq (500u/ μ l) can perform consistently well on a broad range of templates (including both GC and AT rich). PCR BIO Taq DNA Polymerase has 5'-3' exonuclease activities, but no 3'-5' exonuclease (proofreading) activity. The enzyme has the same error rate as wild-type taq DNA polymerase, approximately 1 error per 2.0×10^5 nucleotides incorporated. PCR products generated with PCR BIO Taq DNA Polymerase are A-tailed and may be cloned into TA cloning vectors. PCR BIO Taq (500u/ μ l) is a lyophilisation ready polymerase for incorporation into diagnostic kits and point of care diagnostic devices.

Component	2000 units	20000 units	200000 units
PCR BIO Taq (500u/ μ l)	1x 4 μ l	1x 40 μ l	1x 400 μ l
HS Additive (4 μ g/ μ l)	100 μ l	1mL	10mL

Shipping and Storage

On arrival the kit should be stored at -20°C. Avoid prolonged exposure to light. If stored correctly the kit will retain full activity for 12 months. The kit can be stored at 4°C for 1 month. The kit can go through 30 freeze/thaw cycles with no loss of activity.

Limitations of product use

The product may be used only for in vitro research purposes.

Technical support

For technical support and troubleshooting please email technical@pcrbio.com the following information:

- Amplicon size
- Reaction setup
- Cycling conditions
- Screen grabs of gel images

www.pcrbio.com

Important considerations

Enzyme Concentration: For fast PCR we recommend using between 1 and 5 units of PCRBIQ Taq DNA polymerase per 50µl reaction.

Hotstart Additive: For optimal hotstart the ratio of HS additive to PCRBIQ Taq must be 200ng HS additive/Taq unit. The hotstart additive is supplied in 10mM Tris-HCl, pH8.0. The additive is ready to lyophilise, it does not contain glycerol.

Template: For eukaryotic DNA use between 5ng and 500ng per reaction, for cDNA use below 100ng per reaction.

Primers: Primers should have a predicted melting temperature of around 60°C, using default Primer 3 settings (<http://frodo.wi.mit.edu/primer3/>). The final primer concentration in the reaction should be between 0.2µM and 0.6µM.

Annealing: We recommend performing a temperature gradient to experimentally determine the optimal annealing temperature. Alternatively, we recommend a 55°C annealing temperature then increase in 2°C increments if non-specific products are present.

Extension: Optimal extension is achieved at 72°C. The optimal extension time is dependent on amplicon length and complexity of template. 15 seconds per kilobase (kb) is recommended for amplification from eukaryotic DNA for amplicons between 1kb and 6kb. For shorter amplicons a 1 second extension is sufficient.

10x Reaction buffer: We suggest preparing 10x reaction buffer (100mM Tris-HCl pH8.3 @ 25 °C, 500mM KCl, 30mM MgCl₂, 10mM dNTPs, 1.5% preferred detergent). Filter sterilize.

Enhancers: Standard PCR enhancers may be added to the reaction.

Reaction setup

1. Prepare a master mix based on the following table:

Reagent	50µl reaction	Final concentration	Notes
10x Reaction buffer	5.0µl	1x	Prepared as above
Forward primer (10µM)	2.0µl	400nM	See above for optimal primer design
Reverse primer (10µM)	2.0µl	400nM	See above for optimal primer design
Template DNA	<100ng cDNA, <500ng genomic	variable	See above for template considerations
PCRBIQ Taq DNA polymerase (500u/µl)	0.002µl to 0.01µl	1 to 5u per reaction	
HS Additive (4µg/µl)	0.05µl to 0.25µl	Use 200ng per unit	Correct ratio of HS additive:Taq is essential for optimal hotstart
PCR grade dH ₂ O	Up to 50µl final volume		

2. Cycle using conditions based on the following table:

Cycles	Temperature	Time	Notes
1	95°C	1min	Initial denaturation
40	95°C	15 seconds	Denaturation
	55°C to 65°C	15 seconds	Anneal
	72°C	1 to 90 seconds	Extension (15 seconds per kb)