

EVOLUTION

MORE THAN JUST BASIC SEMEN ANALYSIS



MICROPTIC
AUTOMATIC DIAGNOSTIC SYSTEMS


MICROPTIC HAS SPECIALIZED IN THE ANALYSIS OF SEMEN FOR OVER 20 YEARS




Microptic are the innovators of the world's most advanced CASA (Computer Aided Semen Analysis) system. With dynamic and highly qualified staff, we continue to develop our products in collaboration with the world's leading research centres and laboratories.

SCA[®], THE MOST ADVANCED AND INTEGRATED MODULAR CASA

Automatic measuring system for the analysis of semen samples, following WHO¹ criteria or others.




MOTILITY AND CONCENTRATION




MORPHOLOGY



DNA FRAGMENTATION



VITALITY



ACROSOME REACTION

WHY USE A CASA SYSTEM?

RELIABLE AND DEMONSTRABLE DIAGNOSTICS

Many scientific papers have revealed the homogeneous results obtained with CASA², in contrast to the often unreliable outcome from conventional analysis³.

IMPROVE THE QUALITY OF YOUR LABORATORY PROCEDURES

SCA[®] will assist in the implementation of the objective analysis and quality controls in your laboratory, in addition it produces all of the information needed for the norm ISO 15189:2013 accreditation.

STANDARDIZATION AND TRACEABILITY

The SCA[®] greatly improves the ability of the laboratory to compete through the consistency of the analysis by recording image sequences and results. Measurement errors due to human factors are minimized in the analysis of human sperm.

WIDELY VALIDATED

Various studies have shown that the SCA[®] provides reliable, linear and accurate results with less variability than manual methods, giving higher predictive values when diagnosing fertility problems⁴.

¹ (WHO laboratory manual for the examination and processing of human semen Fifth edition 2010)

² (Maree L., du Plessis S.S., Menkveld R., van der Horst G, Human Reproduction 2010)

³ (Mukhopadhyay D, Varghese AC, Nandi P, Banerjee SK, Bhattacharyya AK. Andrologia, 2010)

⁴ (Dearing C G, Kilburn S, Lindsay K. S., Human Fertility, 2013)

ADAPTATION AND FLEXIBILITY



We have numerous SCA[®] editions with optional modules to help you to find the optimum configuration for your laboratory.

THE MOST COMPLETE CASA SYSTEM



SCA[®] EDITIONS

Several editions with options adapted to give users a choice that will suite their needs:

Lite

Reduced version, just with the basic results for small labs.

Human

Full version, with all the features included for human andrology labs.

Research

Full version to analyse human and animal sperm focused to research centres.

ADDITIONAL MODULES

Enable to complete your SCA[®] CASA system:

SCA[®] dataShare

It enables the internal SCA[®] database sharing and the connection with any LIS system.

SCA[®] Stage Controller

Fully automation using motorized stage.

SCA[®] Manual counter

Manual counter of any biological sample.

COMPLEMENTARY SYSTEMS

Unitary systems that permit to capture, view or edit SCA[®] files:

SCA[®] Capture

Image capture software to use in connection with a main SCA[®] station.

SCA[®] Editor



Software that allows the opening of SCA[®] sessions enabling to edit or reanalyze.

SCA[®] Viewer

Free software that allows to display SCA[®] analysis in any computer.

SCA[®] EVOLUTION: A STATE-OF-THE-ART CASA

Since its launch in 1997, the SCA[®] sperm analyser has been in continuous development, improving with each version the detection and analysis algorithms. Combined with the latest state-of-the-art technology **SCA[®]** becomes the ultimate CASA research tool



VITALITY IN BRIGHTFIELD

We developed BrightVit, a new kit to stain in one step to assess quickly the vitality in brightfield (possible to see the HOS). SCA[®] Vitality can automatically analyse in brightfield or fluorescence.

SPERM FUNCTIONAL TEST

Automatic analysis of hyperactive and mucus penetration, information about sperms capacity to fertilize the egg.

DIAGNOSTIC AND TREATMENT

Automatically provides a diagnosis and recommends the best assisted reproduction treatment to follow.

CUSTOMIZABLE DYNAMIC TABLES

Enable the operator to amalgamate any of the database results into subgroups and create customised statistics.

ADAPTED TO THE NEW TECHNOLOGIES

Working with touch screens, multi-screens, leap motion¹ and portable devices.

ACROSOME REACTION ASSAY

Automatic analysis of intact and acrosome reacted sperm with FluoAcro kit.

ADVANCED ANALYSIS

SCA[®] offers several innovations in automatic analysis such as: an intelligent filter for sperm detection; fluorescence analysis for sperm concentration, motility and DNA fragmentation; very low concentration samples analysis; and vacuoles and tails detection in morphology.

AUTOMATION

With a motorized stage, the analysis process is fully automatic enabling the examination of 4 slides, without user intervention.

RESEARCH EDITION


Analysis of several species, export images, videos and measurement data. Designed for centres involved in research that requires a flexible and exhaustive system.

NEW INTERFACE

Fully customizable.

¹www.leapmotion.com

OUR OBJECTIVE: WARRANTY AND TRUST IN THE DIAGNOSTICS




Reference: 1001


Date: 20/12/2011

Doctor: Dr. Johan Malkovic

Code Patient: 001

Patient: Surname, First Name

Reference: 

Code Patient: 

Collection date: 20/12/2011 08:31	Collection method: Laboratory
Volume (ml): 2.50 (Ref: > 2.5 ml)	pH: 7.50 (Ref: > 7.2)
Temperature (°C): 37	Abstinence days: 3 (Ref: 2-5)
Collection difficulties: No	Complete Sample: Complete
Liquefaction: 60min (Ref: > 60min)	Aggregation: Normal
Viscosity: Normal	Agglutination: 1A
Appearance: Normal	Cellular elements: Normal
Treatment:	

Diagnosis:

Asthenozoospermia

Concentration: (20/12/2011 9:32:00)

28.04 M/ml (Ref: > 15 M/ml) 70.11 M/Ejac (Ref: > 39 M/Ejac) 2.50 ml (Ref: > 1.5 ml) Dilution 1:0

WHOS

	Total	%	Conc/Ejac (M)	Ref
Progressive (PR) (A+B)		23.77	16.67	(> 32%)
Non-progressive (NP) (C)	52	19.62	13.76	
Immotile (IM) (D)	150	56.60	39.68	

	Total	%	Conc/Ejac (M)	Ref
Mobile (A+B+C)	115	43.40	30.42	(> 40%)


Velocity

	Total	%	Conc/Ejac (M)
Rapid	57	21.51	15.08
Medium	32	12.08	8.47
Slow	26	9.81	6.88
Immotile (IM)	150	56.60	39.68

WHO4

	Total	%	Conc/Ejac (M)	Ref
Rapid progressive (type A)	39	14.72	10.32	(> 25%)
Medium Progressive (type B)	24	9.06	6.35	
Non progressive (type C)	52	19.62	13.76	
Immotile (type D)	150	56.60	39.68	

1 / 3



CUSTOMIZED REPORTS SERVICE

More than 40 different reports, including Levy-Jennings report for quality control.

ANALYSIS RESULTS ON-LINE

Patients can have access to their results and analysis videos on: www.semenanalysis.net

mySCA

Virtual site for customers with unlimited space to upload sessions, reports and videos. Also permits users to access webinars, tutorials and on-line teaching courses as well as perform an External Quality Control (EQA).

TRAINING IN YOUR LABORATORY

Installation, training and accreditation by specialist.

IMAGES AND GRAPHICS

The user can select the graphics /images that will appear in the printed report.

RESULTS VALIDATION

A specialist can validate all the assessments, automatically send the results to the central database (e.g.: HL7), and print the custom reports.

LABORATORY INFORMATION SYSTEM (LIS) CONNECTIVITY

SCA® is ready to connect with existing databases.

SUPPORT

Technical advice and service, giving the latest program version.

SCA® MODULES AND SYSTEMS:

MAIN MODULES	
SCA® Motility and concentration	Analysis of the motility, concentration and kinematic parameters
SCA® Morphology	Analysis of the morphology and morphometry in stained sperm samples
SCA® DNA Fragmentation	Analysis of the DNA fragmentation with the chromatine dispersion test
SCA® Vitality	Analysis of the vitality under brightfield or fluorescence
SCA® Acrosome reaction	Analysis of intact and acrosome reacted under fluorescence
SCA® Sample Management	Management of the database and reports generator
ADDITIONAL MODULES	
SCA® DataShare	It enables the internal SCA® database sharing and the connection with any LIS system
SCA® Stage Controller	Fully automation using motorized stage
SCA® Manual Counter	Manual counter of any biological sample
COMPLEMENTARY SYSTEMS (To be used in combination with a main analysis system)	
SCA® Capture	Unitary module for image capture
SCA® Editor	Module for analysis modification
SCA® Viewer	Free software for image visualization

MINIMUM REQUIREMENTS:

	SCA® Motility and Concentration		SCA® Morphology	SCA® DNA Fragmentation		SCA® Vitality		SCA® Acrosome reaction
COMPUTER	Desktop or laptop: Operating system Windows 7 SP1, Windows 8 or Windows 10 (32 or 64 bits); Processor: Intel Core i5 or superior; RAM: 4 GB or superior; USB 3.0 port							
CAMERA	Basler Ace acA1300-200uc / Basler Ace acA780-75gc							
MICROSCOPE	Nikon, Olympus, Zeiss or Leica; Trinocular tube 1x, C-mount 1x; Turret condenser and centering telescope							
OBSERVATION METHOD	Positive phase contrast	Fluorescence	Brightfield	Brightfield	Fluorescence	Brightfield	Fluorescence	Fluorescence
OBJECTIVE	10x Ph+		60x (acA1300) 100x oil (acA780)	10x (acA1300) 20x (acA780)		20x		40x Plan Fluor
FILTERS	Green	Long bandpass (EX 330-380; EM 420; DM 400)	Blue	Green	Bandpass (EX 510-560; EM 590; DM 565)	-	Long bandpass (EX 330-380; EM 420; DM 440)	Dual band (EX 387/480; EM 433/530; DM 403/502)
DISPOSABLES, KITS, STAININGS	Disposable counting chambers or Makler		Sperm Stain, SpermBlue, Diff-Quik, Prestained slides GoldCyto-SB, Papanicolau, Shorr	Chromatin dispersion test kits (SCD)	BrightVit	FluoVit	FluoAcro	
	-	FluoCount						



Health
Canada

DISTRIBUTOR:



MICROPTIC
AUTOMATIC DIAGNOSTIC SYSTEMS



Av. Josep Tarradellas, 8, 1º 6ª · 08029 · Barcelona (Spain) · T. +34 93 419 29 10 · F. +34 93 419 94 26
www.micropticsl.com · micropticsl@micropticsl.com