



## TEST PROCEDURE

Fast, easy and highly specific DNA hybridization tests.



## INNO-LiPA™ TEST PRINCIPLE AND MAIN STEPS

1. Denaturation of amplified biotinylated DNA
2. Hybridization with specific oligonucleotide probes immobilized as parallel lines on membrane-based strips
3. Remove non-specific and unbound DNA
4. Incubation with conjugate and substrate resulting in a purple precipitate



## HBV LiPA COMPATIBILITY

- All LiPA HBV assays
  - use the same program on AutoBlot 3000H or Auto-LiPA™ 48
  - have the same manual protocol
  - have identical LiPA reagents
- Amplicon of INNO-LiPA™ HBV Multi-DR can also be used for INNO-LiPA™ HBV Genotyping

## LiRAST™ FOR LiPA HBV:

### EASY AND OBJECTIVE INTERPRETATION SOFTWARE FOR THE DEVELOPED STRIPS

- Scanning mode with integrated calibration
- Choice between data entry model: scanned or manual entry
- Customized reports: standard or summary report, adjustable to your needs
- Saves electronic image of each strip
- User-friendly and customizable interfaces
- Filter management: e.g. 100 samples with genotype A from January till March
- Patient follow-up\*: overview of test performed per patient and follow-up of drug resistance patterns over time

Interpretation									
Result									
Codon	80	173	180/181	204	236	184	194	202	250
Class	WT	WT	MT/WT	MT	WT	WT	WT	WT	MT
Test Info									
Strips									
LiPA HBV Multi-DR v2 Strip:									
Lot #: not supplied / Assay Date: 11/05/2012 on Sheet: 20120511_001 (Position 17).									
Line Reactivity: 1, 2, 3, 8, 15, 20, 29									
LiPA HBV Multi-DR v3 Strip:									
Lot #: not supplied / Assay Date: 11/05/2012 on Sheet: 20120511_001 (Position 18).									
Line Reactivity: 1, 2, 3, 11, 13, 22									

\* Only for CE version

## ORDERING INFORMATION

### PRODUCT

#### Assays

Assays	CE	RUO
INNO-LiPA™ HBV Genotyping	CE	
INNO-LiPA™ HBV Multi-DR	CE	
INNO-LiPA™ HBV Genotyping		RUO
INNO-LiPA™ HBV Multi-DR		RUO
INNO-LiPA™ HBV PreCore		RUO

#### Software

Software	CE	RUO
LiRAST™ for LiPA™ HBV v1	CE	
LiRAST™ for LiPA™ HBV v1		RUO

#### Automation

Automation	CE	RUO
Auto-LiPA™ 48	CE	
AutoBlot 3000H	CE	

### REFERENCES

1. Yang et al. Clinica Chimica Acta 2010;411(23-24):1951-1956
2. Lok et al. J Clin Microbiol 2002;40:3729-3734
3. Chen et al. J Med Virol 2004;74:536-42
4. European Association for the Study of the Liver. J Hepatol 2012;57(1):167-185
5. Janssen et al. Lancet 2005;365:123-129
6. Zöllner et al. Hepatology 2004;39(1):42-50
7. Tuailon et al. PLoS ONE 2012;7(3):e32143
8. Lin CL, Kao JH. J Gastroenterol Hepatol 2011;26(Suppl 1):123-130
9. Martinot-Peignoux et al. Poster presented at EASL, 18-22 April 2012
10. Congly et al. Poster presented at AASLD, 9-13 November 2012
11. Martinot-Peignoux et al. Poster presented at AASLD, 4-8 November 2011
12. Moucari et al. Antivir Ther 2009;14:1183-1188
13. Gish et al. J Viral Hepat 2010;17(1):16-22
14. Jaroszewicz et al. J Hepatol 2010;52(4):514-522
15. Libbrecht et al. J.Clin.Microbiol 2007;45(12):3935-3941
16. Lampertico et al. Hepatology 2005;42:1414-1419
17. Cornberg et al. Z Gastroenterol 2007;45:1-50
18. Lok AS, McMahon BJ. Hepatology 2009;Vol.50.No.3
19. European Association for the Study of the Liver. J Hepatol 2009;50(2):227-42
20. Zoulim et al. J Viral Hepat 2007;14 Suppl 1:29-36
21. Sablon et al. Expert Rev Mol Diagn 2003;3(5):535-47
22. Hussain et al. J Clin Microbiol 2003;41(8):3699-3705
23. Ozasa et al. Hepatology 2006;44(2):326-34
24. Malik et al. PLoS One 2012;7(6):e39028
25. Hayashi et al. Intervirology 2009;52(1):22-8
26. Xiao et al. J Med Virol 2011;83(9):1544-50
27. Ren et al. J Viral Hepat 2010;17(12):887-95
28. Sonneveld et al. Hepatology 2012;56(1):67-75
29. Tacke F, Shirvani-Dastgerdi E. Hepat Mon 2012;12(6):357-60

Fujirebio Europe N.V.  
Technologiepark 6- B-9052 Gent - Belgium  
T +32 (0)9 329 13 29 - F +32 (0)9 329 17 75  
www.fujirebio-europe.com

### ARTICLE NO.

Assays	CE	RUO	ARTICLE NO.
INNO-LiPA™ HBV Genotyping	CE		80691
INNO-LiPA™ HBV Multi-DR	CE		81383
INNO-LiPA™ HBV Genotyping		RUO	80070
INNO-LiPA™ HBV Multi-DR		RUO	81382
INNO-LiPA™ HBV PreCore		RUO	80883


Software	CE	RUO	ARTICLE NO.
LiRAST™ for LiPA™ HBV v1	CE		81387
LiRAST™ for LiPA™ HBV v1		RUO	81386

Automation	CE	RUO	ARTICLE NO.
Auto-LiPA™ 48	CE		80628
AutoBlot 3000H	CE		81149

© March 2014, Fujirebio Europe

The complete portfolio for the molecular identification of the Hepatitis B Virus

INNO-LiPA™ HBV Genotyping INNO-LiPA™ HBV Multi-DR INNO-LiPA™ HBV PreCore



- Fast, accurate and informative
- Easy to perform
- Highly sensitive
- Full patient profiling with 1 method
- Sensitive detection of (mixed) genotypes
- Earliest detection of drug resistance
- Accurate determination of the true HBeAg status
- HBV treatment management

# INNO-LiPA™ HBV Genotyping

Line probe assay for the identification of hepatitis B virus genotypes A to H

## STRIP LAY-OUT



## FEATURES AND BENEFITS

- Clear identification of mixed infections: LiPA detects up to 16.3 % more mixtures than direct sequencing, confirmed by clonal analysis <sup>[1, 2]</sup>
- Sensitive detection of genotypes: LiPA is more sensitive than PCR-RFLP (98.8 % vs 65.0 %) <sup>[3]</sup> and more sensitive than sequencing <sup>[2]</sup>
- Amplicon of INNO-LiPA™ HBV Multi-DR can also be used for INNO-LiPA™ HBV Genotyping


## CLINICAL RELEVANCE OF HBV GENOTYPING

- Genotyping linked with treatment response and part of clinical guidelines <sup>[4]</sup>: HBV genotypes A and B are associated with a better response to pegylated IFN- $\alpha$  than genotypes C and D. <sup>[5]</sup> The mutational pattern of lamivudine drug resistance differs between genotypes A and D. <sup>[6]</sup>
- The correlation between levels of Hepatitis B surface antigen (HBsAg) and the HBV genotype makes both markers valuable tools for prognosis and management of chronic hepatitis B: HBsAg levels differ according to HBV genotype prior to treatment <sup>[7-12]</sup>, during treatment with pegylated IFN- $\alpha$  <sup>[8, 12]</sup> and nucleos(t)ide analogues <sup>[13]</sup> and during the natural course of HBV infection. <sup>[14]</sup>

# INNO-LiPA™ HBV Multi-DR

Line probe assay for the simultaneous detection of hepatitis B virus wild-type motif and resistance associated mutations

## STRIP LAY-OUT



## FEATURES AND BENEFITS

- Relevant, accurate and up-to-date detection of drug resistance associated mutations as well as known compensatory mutations
- LiPA is highly sensitive and specific and offers the earliest timepoints for detection of mutants:
  - 2,8 times greater chance for LiPA to detect a given mutation than direct sequencing at any moment in time <sup>[15]</sup>
  - early detection of genotypic resistance allows rapid and effective HBV suppression <sup>[16]</sup>
- HBV drug resistance testing: part of clinical guidelines <sup>[17-19]</sup>


## CLINICAL RELEVANCE OF HBV DRUG RESISTANCE TESTING

- Antiviral drug resistance in patients with chronic hepatitis B leads to antiviral treatment failure. <sup>[20]</sup>
- Resistance should be identified as early as possible before viral load rebound and before rise of ALT (alanine aminotransferase) levels. <sup>[21]</sup>
- As soon as drug resistance is identified, an appropriate therapy can be initiated with the most effective antivirals now available. <sup>[16]</sup>

# INNO-LiPA™ HBV PreCore

Line probe assay for the simultaneous detection of HBV wild-type motif and mutations in the basal core promotor (BCP) and PreCore regions

## STRIP LAY-OUT



## FEATURES AND BENEFITS

- More sensitive than direct sequencing for the detection of mixed sequences in precore (PC) and basal core promotor (BCP) region <sup>[22]</sup>
- Prediction of increased risk for hepatocellular carcinoma and progression of liver disease <sup>[23-27]</sup>
- Determine cause of HBeAg negativity: treatment response or PC/BCP mutations? <sup>[28-29]</sup>

## RELEVANCE OF HBV PRECORE IN CLINICAL PRACTICE

- Indication on treatment response and duration:
  - Pegylated IFN- $\alpha$  therapy response: patients without PC/BCP mutations are good candidates for interferon treatment, independent of the HBV genotype <sup>[28]</sup>
  - Nucleos(t)ide analogues treatment + PC/BCP mutations + drug resistance mutations might lead to a higher or lower replication fitness of the resistant virus <sup>[29]</sup>
- Risk of fulminant hepatitis and hepatocellular carcinoma: chronic hepatitis B patients infected with PC/BCP mutant viruses are more susceptible to severe hepatitis and acute-on-chronic liver failure (ACLF) <sup>[23-27]</sup>