

# MVE CryoCube™



The MVE CryoCube™ provides a cryogenic shipping option for your biological samples. In keeping with MVE tradition, the CryoCube incorporates creative engineering and simplistic design to allow shipping in any orientation without sacrificing temperature or hold time.

## Features include:

- Charges in under 1 hour
- Operates as a dry shipper
- -150°C temperature
- 5-day holding time
- Weighs less than 9 lbs. charged
- Safe shipping in any orientation
- Safe for international shipping

## CryoCube™

### Maximum Storage Capacity

Number of Canisters	1
Number of 1/2 cc Straws 1 Level Bulk	88
Number of 1/4 cc Straws 1 Level Bulk	182
Number of 1.2 & 2.0 ml Vials bulk	6

### Performance

LN2 Capacity(min) L	2
Static Holding Time days	5

### Unit Dimensions

Neck Opening in. (mm)	1.4 (35)
Overall Height in. (mm)	12 (305)
Outer Width in. (mm)	12 (305)
Inner Canister Height in. (mm)	5.25 (127)
Inner Canister Diameter in. (mm)	1 (35)
Weight Empty lb. (kg)	5.25(1.8)
Weight Charged lb. (kg)	9 (4)

90 Day Standard Warranty

Go to [www.chartbiomed.com](http://www.chartbiomed.com) for more information.

Chart Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

