

CryoPur™

Cryopreservation Solutions

CryoPur™: 100% DMSO

Sterile filtered DMSO is available in a variety of packages, concentrations, and volumes designed for lab use. 100% and buffered DMSO solutions are packaged in polymeric syringes or vials; both packages are water-clear, and shatter-resistant. Our DMSO products are CE-marked and certified to USP and Ph. Eur. standards

CryoPur™ vials with white caps contain 100% DMSO, and are available in 10, 50, 70 and 100ml fill volumes. Caps are flip off/tear off seal.

CryoPur-D™ DMSO/Dextran Solutions

Sterile, filtered DMSO/Dextran solutions are available in a variety of concentrations and formulations. CryoPur-D™ vials with green flip/tear off caps contain 55% w/v DMSO, 5% with Dextran-40, remainder water for injections.

Available in Syringes and Vials

DMSO and DMSO/Dextran filled sterile syringes are fitted with a stopcock and are shipped with the actuator blocked for safety. To use merely connect to your system. They are available in several configurations for easy, aseptic transfer.



CryoPur™ Cryopreservation Solutions

Ordering Information

CryoPur™: 100% DMSO



Product Description	Stock Item	Stock number	European Ref. Number
Vial, 10ml, white 'flip/tear up' cap	YES	CP10	V10010
Vial, 50ml, white 'flip/tear up' cap	YES	CP50	V10050
Vial, 70ml, white 'flip/tear up' cap	YES	CP70	V10070
Vial, 100ml, white 'flip/tear up' cap	YES	CP100	V10100
Syringe, 10ml, with stopcock and 15cm SCD tubing		SP10	S10010
Syringe, 50ml, with stopcock and 15cm SCD tubing		SP50	S10050

CryoPur-D™: 55% w/v DMSO, Buffered with 5% w/v Dextran-40

Product Description	Stock Item	Stock number	European Ref. Number
Vial, 7.2ml, green 'flip/tear up' cap	YES	CD 07	V504007
Polymer Syringe, 10ml, with stopcock and 15cm SCD tubing		SD 10	S504010

Storage: 100% DMSO: +18 to 25°C; DMSO/Dextran: +2 to 8°C

Packaging: 10ml vials 12 to a case, larger vials 6 to a case, syringes sold individually

Custom Sizes: Custom packages and formulations possible, please enquire

USA: For Research use only. 8-12 week lead times on custom items.

OriGen

BIOMEDICAL

Americas & Asia

7000 Burlleson Rd, Bldg D
Austin, TX USA 78744-3202
Tel: +1 (512) 474 7278
Fax: Toll Free (888) 812 8411 or
+1 (512) 617 1503

Europe & Middle East

Landskronavägen 25A
252 32 Helsingborg, Sweden
Tel. + 46 46 772 433
Fax. +46 44 785 1099

www.origen.com