

Incubation Trays

Convenient, easy to use trays for staining and washing gels and membranes

These specially designed **incubation trays** are ideal for staining and washing electrophoresis gels and membranes.


Choose from three designs, all available in a range of sizes:

- **Opaque** – for light-sensitive applications
- **Transparent** – for easy monitoring of colorimetric staining
- **Traditional** – perfect for blot washes and incubations




Advantages

- **SMOOTH INTERIOR** – protects membranes and gels from scratches
- **ATTACHED LID** – to protect experiments from dust or debris that can cause speckles on blot images
- **CONVENIENT** – multiple sizes available for different sizes of gels and blots
- **SAVE** – minimize antibody and buffer usage by using appropriately sized trays



Easy visualization of staining



Tray lids protect from dust and debris that cause speckles on Western blots



Save antibody and buffer by using the smallest tray that can hold your blot



Beautiful Western blots every time

Advansta Solutions

Staining with a light-sensitive dye or doing a fluorescent Western?

Use Advansta's opaque incubation trays

- Protect your blot or gel from light
- No need to cover your tray with foil
- No need to keep your shaker in a darkened room
- Use for hybridization and washes of fluorescent microscope slides

Not sure when staining or destaining is complete?

Use Advansta's transparent incubation trays

- Easily observe staining
- Photograph your gel right in the tray

Incubation Trays

Convenient, easy to use trays for staining and washing gels and membranes



Ordering Information

Catalog Number	Product	Size	Size
L-07034-005	Incubation Tray, Transparent body and lid	5x7 cm	5/pack
L-07038-005	Incubation Tray, Opaque black body and lid	5x7 cm	5/pack
L-07031-005	Incubation Tray, Traditional red body with clear lid	5x7 cm	5/pack
L-07035-004	Incubation Tray, Transparent body and lid	6x9 cm	4/pack
L-07039-004	Incubation Tray, Opaque black body and lid	6x9 cm	4/pack
L-07032-005	Incubation Tray, Traditional black body with clear lid	6x9 cm	5/pack
L-07036-003	Incubation Tray, Transparent body and lid	9x11 cm	3/pack
L-07040-003	Incubation Tray, Opaque black body and lid	9x11 cm	3/pack
L-07033-005	Incubation Tray, Traditional blue body with clear lid	9x11 cm	5/pack
L-07037-002	Incubation Tray, Transparent body and lid	10.5x15.5 cm	2/pack
L-07041-002	Incubation Tray, Opaque black body and lid	10.5x15.5 cm	2/pack

Related Products

K-12045-C20	WesternBright™ ECL Western Blotting HRP Substrate	Trial kit size	20 ml
K-12045-D20	WesternBright™ ECL Western Blotting HRP Substrate	For 2000 cm ² membrane	200 ml
K-12045-D50	WesternBright™ ECL Western Blotting HRP Substrate	For 5000 cm ² membrane	500 ml
R-05071-500	Goat-anti-mouse HRP-conjugated secondary antibody		500 µl
R-05072-500	Goat-anti-rabbit HRP-conjugated secondary antibody		500 µl
L-08001-010	Low-Fluorescence PVDF Transfer Membrane	7x9 cm	10 sheets
L-08002-010	Nitrocellulose Transfer Membrane 0.45 µm	8x10 cm	10 sheets
L-08003-010	Nitrocellulose Transfer Membrane 0.22 µm	8x10 cm	10 sheets

Advansta Corporation

1505 Adams Drive, Suite B1 | Menlo Park, CA 94025

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com

Product information: www.advansta.com/Incubation_trays

Copyright © 2010–2014 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. WesternBright™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

www.advansta.com

