

# Pre-cut membranes for immunoblotting

## Convenient Nitrocellulose and PVDF membranes for Western blotting

Advansta offers convenient **pre-cut membranes** for Western blotting. Sized for protein minigels, these membranes are available in both nitrocellulose and low-fluorescence PVDF. Save time and avoid concerns about membrane contamination from dirty scissors or gloves with quality-controlled, ready-to-use membranes.



### Advantages

- **PERFORMANCE** – quality controlled for background-free Western blots
- **CONVENIENCE** – save time, and avoid concerns about membrane contamination
- **FLEXIBILITY** – choose the right membrane for your Western blotting application

### Compatible with gels 10 x 10 cm and smaller, including

Manufacturer	Product Name
Bio-Rad	Mini-PROTEAN® Tetra system, Mini-PROTEAN and Ready Gel® precast gels
GE Healthcare	miniVE system, Hoefer™ SE 250 and SE 260 Mighty Small II systems
Invitrogen	NuPAGE® Novex® and Novex ZOOM® precast gels, XCell6™ MultiGel system, XCell Surelock™ system
Pierce	Precise™ Protein Gels
Thermo Electron	EC120 Mini Vertical Gel System

### Ordering Information

Catalog Number	Product		Size
L-08001-010	Low-Fluorescence PVDF Transfer Membrane	7 x 9 cm	10 sheets
L-08012-010	Low-Fluorescence PVDF Transfer Membrane	10 x 15 cm	10 sheets
L-08014-010	Low-Fluorescence PVDF Transfer Membrane	13 x 18 cm	10 sheets
L-08002-010	Nitrocellulose Transfer Membrane 0.45 µm	8 x 10 cm	10 sheets
L-08003-010	Nitrocellulose Transfer Membrane 0.22 µm	8 x 10 cm	10 sheets

### Related Products

K-12021-010	WesternBright™ MCF fluorescent Western blotting kit (Contains: secondary antibodies, washing solution, blocking solution, pre-cut PVDF membranes and background quenching sheets)		10 blots
K-12042-D10	WesternBright™ Quantum Western Blotting HRP Substrate	For 1000 cm <sup>2</sup> membrane	100 ml
K-12045-D20	WesternBright™ ECL Western Blotting HRP Substrate	For 2000 cm <sup>2</sup> membrane	200 ml

### Advansta Corporation

1505 Adams Drive, Suite B1 | Menlo Park, CA 94025

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: [sales@advansta.com](mailto:sales@advansta.com)

Product information: [www.advansta.com/Transfer\\_membranes](http://www.advansta.com/Transfer_membranes)

Copyright © 2010-2014 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. WesternBright™ and LucentBlue™ are trademarks of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

[www.advansta.com](http://www.advansta.com)

