

MODEL	S.C. 2/1V	S.C. 4/2V	S.C. 4/3V	S.C. 20/12V	XC 20/3V	Mini Moover	Cryo	Cryo	Cryo	IATA
							Moover	Shipper	Shipper XC	
Pojemność [l]	1.5	3.6	4.3	20.5	6.2	2.9	4.2	8.5	10.0	10.0
Szybkość parowania azotu l/24h	0.19	0.26	0.20	0.09	0.35	0.20	0.35	0.85	0.70	0.70
Statyczny czas pracy [ilość dni]	8	13	21	60	16	14	12	10	14	14
SPECYFIKACJA										
Szerokość otworu [mm]	35	70	51	55	96.7	35	97	216	216	216
Wysokość całkowita [mm]	343	468	492	652	559	495	558	546	584	584
Średnica całkowita [mm]	184	222	222	368	368	184	464	369	381	381
Wysokość kanistra [mm]	127	279	279	279	279	279	279	279	317.5	222.25
Średnica kanistra [mm]	31	67	46	38	80	31	79	142 x 167	216	187.325
Waga netto [kg]	2.7	4.5	5.3	12.8	11.3	4.1	13.8	11.8	15.0	15.0
Waga brutto [kg]	4.1	9.1	9.5	29.3	27.3	7.0	39.5	24.7	33.1	33.1
Ilość kanistrów	1	1	1	6	4+1 Canister	1	7	1Rack	-	Secondary Container
<b>Ilość słomek ½, 10/gobl.</b>	-	280	120	540	2500/2000	60	3080	-	-	-
Ilość słomek ½, luzem	88	44	210	780	3750/3000	88	4354	-	-	-
Ilość probówek 1.2&2ml, 5/gobl.	-	95	40	150	675/560	20	945	-	-	-
Ilość probówek 1.2&2ml, 6/gobl.	-	102	48	180	840/672	24	1134	500	966 (luzem)	-
Ilość worków 4R9953 (250ml)	-	-	-	-	-	-	-	10	10	-