

MODEL	XC 20 Signature	XC 21/6	XC 22/5	XC 32/8	XC 33/22	XC 34/18	XC 34/18 Plus	XC 43/28	XC 47/11-6SQ	XC 47/11-6	XC 47/11-10
Pojemność [l]	20.5	21.0	22.4	32.0	33.4	34.8	34.8	42.2	47.4	47.4	47.4
Statyczna szybkość parowania azotu [l/24h]	0.090	0.35	0.35	0.35	0.14	0.18	0.40	0.14	0.39	0.39	0.39
Normalny czas pracy [ilość dni]	142	38	40	57	154	123	136	193	76	76	76
SPECYFIKACJA											
Szerokość otworu [mm]	55.4	89	97	97	70	89	89	70	127	127	127
Wysokość całkowita [mm]	652	438	559	546	660	675	675	670	673	673	673
Średnica całkowita [mm]	368	464	368	464	464	464	464	508	508	508	508
Wysokość kanistra [mm]	279	127	279	279	279	279	279	279	-	279	279
Średnica kanistra [mm]	41.9	70	79	67	56	71	71	56	-	102	71
Waga netto [kg]	10.0	10.9	11.8	13.6	15.4	15.4	20.8	16.4	19.0	19.0	19.0
Waga brutto [kg]	26.5	27.9	30.0	39.5	42.5	43.5	47.3	50.5	54.6	54.6	54.6
MAKSYMALNA POJEMNOŚĆ											
Ilość kanistrów	6	9	6	9	6	6	6	6	6 sq.	6	10
Ilość słomek ½ 10/gobletka	660	N/A	2,400	2,520	1,260	2,100	2,100	1,260	-	3,900+	3,200
Ilość słomek ½, luzem	879	3,870	3,666	3,960	1,764	3,300	3,300	1,764	-	6,200	5,000
Ilość probówek 1.2&2ml, 5/gobletka	210	N/A	810	855	360	630	630	360	-	1,320	1,050
Ilość raków na 25 probówek	-	-	-	-	-	-	-	-	750	-	-